

Chapter 13

Cell Culture Techniques

After reading this chapter, you should be able to:

- Compare the advantages and disadvantages of using in vitro cell culture techniques with in vivo techniques
- Describe various cell culture preparations used to examine nervous system function in vitro
- Describe techniques for manipulating/influencing cells in culture

Techniques covered:

- **Tools and reagents used in cell culture:** equipment, media
- **Types of cultured cells:** immortalized cell lines, primary cell and tissue culture, stem cells, brain organoids
- **Manipulating/influencing cells in culture:** transfection, coculture, pharmacology, antibody interference

Cell culture is the process by which a scientist grows and maintains cells under carefully controlled conditions outside of a living animal. There are many reasons why this approach is desirable. Examining the nervous system in vitro (“within glass”) allows scientists to simplify the cellular environment, providing greater control over experimental manipulations and reducing potentially confounding interactions with other biological systems. In vitro tools and techniques make experiments possible that would otherwise be difficult or impossible to perform (or interpret) in intact organisms, such as imaging of subcellular structures in living cells or performing reproducible assays with the exact same number of cells. In vitro experiments also tend to be faster, less expensive, and require fewer animals than experiments performed in vivo.

However, cell culture experiments also raise interesting questions: Does a cell in a culture dish truly behave like it would inside the brain or body? What does it mean to culture neurons, cells distinguished by their ability to communicate with other cells in complicated networks, in relatively isolated conditions such that they do not form synaptic connections? Are immortalized cell lines, cells that continue to divide indefinitely, actually good models for cells that degenerate and perish in the brain? It is important to keep these questions in mind when designing in vitro experiments and interpreting results.

This chapter surveys approaches to studying the nervous system using cell culture techniques. We start with a quick summary of the equipment and reagents necessary to maintain cells outside of living organisms. Next, we describe categories of cells used in culture experiments: immortalized cell lines, primary cell cultures, stem cells, and brain organoid systems. Finally, we describe some of the methods scientists can use to manipulate and influence cells in cultured environments.

CELL CULTURE EQUIPMENT AND REAGENTS

Specialized equipment and reagents are necessary to provide cultured cells with an environment that can support their continued growth and health outside a living organism. Most of these supplies are used to artificially mimic the endogenous, *in vivo*, cellular environment. Most equipment and reagents in cell culture labs are used to continually supply cells with oxygen, nutrients, growth factors, and other elements necessary to keep cells alive.

Other tools and reagents prevent contamination. Cells *in vivo* have the benefit of an active immune system to prevent contamination from bacteria, fungi, and other microorganisms. Cells in the brain have the added benefit of the blood–brain barrier to further prevent contamination. In culture conditions, cells are incredibly vulnerable, and scientists must take great care in avoiding contamination, especially because microorganisms are ubiquitous in the environment and can easily penetrate cell culture plates.

Equipment

Although individual laboratories study different cell lines and ask different scientific questions, cell culture rooms tend to contain the same fundamental pieces of equipment. This equipment includes

- *Biosafety Hood.* A biosafety hood or laminar flow cabinet is used to prevent contamination by microorganisms (Fig. 13.1). When not in use, these chambers are often illuminated with UV light to sterilize exposed surfaces. Just before use and throughout experiments, scientists spray all surfaces with 70% ethanol to provide further decontamination. During experiments, scientists always keep reagents and cell culture flasks/plates covered until they are ready for use; when exposed, a scientist should never pass their gloved hands or other equipment over open bottles or, even worse, the cells themselves, as any microorganisms present on gloves or tools could contaminate the media.
- *Cell Incubator.* Cell incubators house and store culture flasks and plates. These incubators maintain an appropriate temperature, humidity, and gas concentration to mimic endogenous conditions. They are usually set at 37°C with 5% CO₂ levels. Although 5% CO₂ is not the level experienced in the body *in vivo*, the carbon dioxide maintains the pH of buffers in the growth media at proper physiological levels.




FIGURE 13.1 Typical cell culture room. A CO₂ incubator is used to keep the cells at an ideal temperature and provide CO₂ for the buffering solutions in the culture media. A biosafety hood is used to prevent contamination from microorganisms during handling.

- *Specially prepared cell culture flasks/plates.* Tissue culture flasks and plates come in many varieties depending on the needs of the investigator. Flasks typically range in sizes of 12.5 cm² to large, multitiered chambers with surface areas of 1800 cm². Multiwelled plates are often used for tissue culture experiments, as these wells can hold an exact number of cells, and it is easy to keep track of various experimental conditions across wells. Plates are typically sold with 6, 12, 24, 48, 96, or 384 wells. Many cells, including neurons, need to adhere to a substrate to grow, so tissue culture dishes often need to be specially prepared to provide a substrate for the cells. Tissue culture dishes or glass coverslips are often coated with attractive amino acids, such as lysine or ornithine, or extracellular matrix components like collagen and laminin. Other materials such as hydrogels or polyacrylamide gels are sometimes used as softer substrates on which cells can grow, and these softer environments are often permissive for cell growth or improved cell health.
- *Water bath.* A water bath is often set at 37°C and is used to quickly warm cell culture media and other reagents stored at 4°C just prior to being added to cell culture flasks/plates. If a scientist does not warm media before adding it to cells, the cells can be shocked by the abrupt decrease in temperature and die.
- *Microscope.* Microscopes are used in most tissue culture rooms for routine observation of cell culture flasks/plates to inspect the health and confluence of cells. **Confluence** refers to the percentage of the surface of the bottom of

the plate covered by cells. Most cell lines should never become 100% confluent, as cells in dense populations tend to inhibit each other's growth, as well as quickly drain nutrients from the culture media.

- *Refrigerator.* A refrigerator (a 4°C incubator) maintains cell culture media and other reagents when not in use. The growth factors and antibiotics in culture media degrade over time, but can last for weeks if stored at 4°C.

Culture Media

Growth media is critical to cell culture experiments, supplying nutrients (amino acids and vitamins) and a source of energy (glucose) for cells. Growth media can vary in pH, nutrient concentration, and the presence of growth factors or other biologically relevant components. To survive, cells must be bathed in an isotonic fluid that has the same concentration of solute molecules as inside the cell. The media is buffered to maintain a compatible pH (usually 7.4, though there are some cell-specific variations). In the 5% CO₂ environment of culture incubators, bicarbonate buffer maintains the physiological pH as well as provides nutritional benefits to the cells. HEPES (4-(2-hydroxyethyl)-1-piperazineethanesulfonic acid) is added to the culture media for extra buffering capacity when the cells in a culture experiment require extended periods of manipulation outside a CO₂ incubator.

Serum is often added to culture media for its ability to promote survival through undefined mixtures of growth factors, hormones, and proteins, like PDGF (platelet-derived growth factor), insulin, and transferrin. However, the use of serum is nonphysiological as most proteins found in serum do not normally cross the blood–brain barrier. Thus, while adding serum to the culture media allows for robust cell survival, it can create cells with properties that are distinct from cells *in vivo*. For example, serum-exposed astrocytes divide rapidly, are fibroblast-like in appearance, and lack the highly branched processes typical of *in vivo* astrocytes (Fig. 13.2).

For stricter control over the cellular environment, investigators use serum-free supplements, such as N2 or B27, that contain known formulations of survival factors. A typical growth medium in neuronal cultures is NeuroBasal, which provides optimized amino acids and nutrients to cultured neurons, supplemented with N2 or B27. Many serum-free glial cultures require specific growth factors to keep the cells alive and nonreactive.

Adding or removing specific ingredients to or from the media can influence cellular behavior. For example, to maintain a neural progenitor pool that continues to divide in culture, scientists add growth factors such as epidermal growth factor (EGF) and fibroblast growth factor (FGF) to prevent differentiation. Specific recipes for culture media are tailored to provide optimal conditions for promoting the health and proper physiological function of the cultured cells and maintaining cellular behavior.




FIGURE 13.2 Serum-treated astrocytes behave differently from quiescent astrocytes. (A) Astrocytes cultured in the absence of serum retain their characteristic morphology and properties. (B) Addition of serum to cultured astrocytes leads to increased proliferation and abnormal morphology. *Reprinted with permission from Foo, L. C. et al. (2011). Development of a method for the purification and culture of rodent astrocytes. Neuron: Cell Press.*

Now that we have described some of the equipment and reagents necessary for cell culture, we survey different categories of cells commonly used in cell culture experiments.

IMMORTALIZED CELL LINES

Immortalized cell lines are either tumorous cells that do not stop dividing or cells that have been artificially manipulated to proliferate indefinitely and can thus be cultured over several generations (Table 13.1). Because immortalized cells continuously divide, they eventually fill up the dish or flask in which they grow. By **passaging** (also known as **splitting**), scientists transfer a fraction of the multiplying cells into new dishes to provide space for continuing proliferation. Common immortalized cell lines include **human embryonic kidney 293T (HEK-293T) cells** and **HeLa cells**.

There are many advantages to using immortalized cell lines. Because these lines are used by many different labs in various experimental contexts, they are well characterized. Furthermore, they are homogeneous, genetically identical populations of cells, allowing for consistent and reproducible results. Immortalized cells tend to be easier to culture than cells used in primary cultures in that they grow more robustly and do not require extraction from a living animal. Also, because they grow quickly and continuously, it is possible to extract large amounts of protein for biochemical assays (Chapter 14). It is also possible to create cell lines that continuously express a gene of interest, such as a fluorescently tagged or mutant version of a protein.

The major disadvantage to using immortalized cells is that these cells cannot be considered “normal,” in that they divide indefinitely and sometimes express unique gene patterns not found in any cell type *in vivo*. Therefore, they

TABLE 13.1 Commonly Used Immortalized Cell Lines.

Cell Line	Origin and Cell Type	Comments
HEK-293T	Human embryonic kidney cell	Easy to transfect and manipulate; commonly used as an expression system to study signaling and recombinant proteins
HeLa	Human epithelial cell	From a cervical cancer in a human patient named Henrietta Lacks; able to grow in suspension (i.e., grow without adhering to the bottom of plates)
COS	Monkey kidney cell	Efficiently transfected; commonly used as an expression system for high-level, short-term expression of proteins
3T3	Mouse embryonic fibroblast	Robust and easy to handle; contact inhibited; stops growing at very high densities
MDCK	Dog kidney epithelial cell	Polarized with distinct apical and basal sides, used in studying trafficking
CHO	Chinese hamster ovary	Useful for stable gene expression and high protein yields for biochemical assays; commonly used as an expression system for studying cell signaling and recombinant proteins
S2	<i>Drosophila</i> macrophage-like cells	Well-characterized <i>Drosophila</i> cell line; highly susceptible to RNAi treatment
PC12	Rat pheochromocytoma chromaffin cell	Neuron-like, derived from a neuroendocrine adrenal tumor; can differentiate into a neuron-like cell in the presence of NGF
Neuro-2a (N2a)	Mouse neuroblastoma	Model system for studying pathways involved in neuronal differentiation; can be driven to differentiate by cannabinoid and serotonin receptor stimulation
SH-SY5Y	Human neuroblastoma, cloned from bone marrow	Grows as clusters of neuroblast-type cells with short, fine neurites; can be dopaminergic, noradrenergic, acetylcholinergic, glutamatergic, adenosinergic

might not have the relevant attributes or functions of typical cells. Also, after several passages, cell characteristics can change and become even more different from those of a normal cell. Thus, it is important to periodically validate the characteristics of cultured cells and not use cells that have been passaged too many times.

Immortalized cell lines of neuronal origin can be used to study properties unique to neurons. Scientists have used neuronal cell lines to investigate processes that occur during differentiation in neurons, such as axon selection, guidance, and growth. However, most neuronal immortalized cell culture models are derived from tumors and are sometimes genomically abnormal. One popular neuronal cell line, called PC12, is a rat pheochromocytoma cell line derived from an adrenal gland tumor. The addition of nerve growth factor causes PC12 cells to reversibly differentiate into a neuronal phenotype. These cells can synthesize dopamine, norepinephrine, and acetylcholine. PC12 cells have been used to study molecular phenomena associated with neuronal differentiation and have even been used in experiments to replace dopaminergic neurons in an animal model of Parkinson's disease. Neuroblastoma cell lines, like mouse Neuro2A, also express neurotransmitters and have been used in electrophysiology and neurodevelopment studies.

As useful as immortalized cell lines of neural origin can be for certain experiments, they show obviously abnormal traits, such as the unusual combination of neurotransmitters they produce (no normal neuron produces dopamine, norepinephrine, and acetylcholine in the same cell!). Therefore, it is advantageous, when possible, to use primary cultured cells—cells extracted from living animals.

PRIMARY CELL AND TISSUE CULTURE

Primary cell culture and **primary tissue culture** are derived from tissue removed directly from a living animal rather than immortalized cells that divide indefinitely. Primary tissue culture allows scientists to directly investigate cells of interest in a carefully controlled in vitro environment. These experiments can be especially advantageous for a comparison between cells derived from a wild-type animal and cells derived from a genetically modified animal.

However, there are some disadvantages to using primary cultures instead of immortalized cell lines ([Table 13.2](#)). Unlike immortalized cell lines, primary cultured cells have a limited lifetime. The age of the animal source influences the health and robustness of the cell culture: tissue from younger, embryonic, or early postnatal animals survives better and tends to be healthier than tissue from older animals. Also, a population of primary cells will always be more heterogeneous than a culture of immortalized cells, no matter how careful the scientist was in extracting and purifying the cells of interest.

There are three main categories of primary tissue culture: slice cultures, explant cultures, and dissociated cultures ([Fig. 13.3](#)).

Slice Cultures

Slice cultures maintain the structure and organization of the brain in a relatively thick (250–400 μm) section of brain tissue cut using a vibratome

TABLE 13.2 Immortalized Cell Lines Versus Primary Cultures.

	Advantages	Disadvantages
Immortalized cell lines	<ul style="list-style-type: none">• Easier to use (grow, passage, transfect, etc.)• Homogeneous population• Well characterized and consistent between labs• Can create stable cell line expressing gene of interest	<ul style="list-style-type: none">• May not have the same properties as neurons or primary cell type of interest• Cell lines do not exist for many neural cell types
Primary cell culture	<ul style="list-style-type: none">• Relevant cell type, physiology, and circuitry• Protocols for culturing many different types of cells from the nervous system	<ul style="list-style-type: none">• Can have heterogeneous populations with high variability• Time-consuming protocols to harvest cells from living organisms

(Chapter 6). Slices provide greater access to and visibility of deep subcortical structures, such as the hippocampus and thalamus, which are difficult to access in vivo (Fig. 13.3). These slices can either be **acute cultures** that are used immediately (over a period of hours) or **organotypic slice cultures** that are maintained over multiple days. Acute slices are typically used for short-term electrophysiology experiments, while organotypic cultures are used to observe structural and morphological changes such as neuronal migration, axon outgrowth, synapse formation, myelin formation, or activation of microglia.

Explant Cultures

Explant cultures consist of intact fragments of tissue (Fig. 13.3). Although they do not necessarily preserve the precise organization and connectivity of the endogenous nervous system, explants contain the same mixture of cell types and some connectivity remains. Explants are often used in coculture assays in which two different cell culture models are present in the same culture dish (see below).

Dissociated Cell Cultures

In **dissociated cultures**, neural tissue is separated into individual cells that are then grown on two-dimensional coated glass coverslips or within three-dimensional substrates (Fig. 13.3). After extracting an animal’s brain, tissue from specific brain regions can be further microdissected (e.g., cortical, hippocampal). This region is then mechanically or enzymatically digested within a liquid suspension to separate individual cells from each other. The neurons




FIGURE 13.3 Primary cell and tissue culture preparations. After a brain is removed from an animal, it can be cut into thin sections (usually using a vibratome) and kept alive in slice culture. Alternatively, either a brain slice or a brain region can be finely dissected to produce an explant culture. An explant culture can be enzymatically digested into single cells to produce a dissociated cell culture.

are then removed from the suspension and plated onto a substrate on which the cells can attach and grow.

Neurons or glial cells dissociated from different regions of the brain retain their initial identities. The morphological, molecular, and physiological properties of cell populations present in culture often correspond closely to the characteristics of the cell population present in the region of origin in a living organism. With the proper growth factors and care, it is possible to maintain a dissociated neuronal or glial culture for weeks, during which time cells acquire properties of mature neurons and glia. Neurons, for example, develop characteristic axons and dendrites, form synapses with one another, and express receptors and ion channels specific to their cell types, even producing spontaneous electrical activity. Dissociated neuronal cell cultures have been used to study neurite outgrowth, synapse formation, and electrophysiological properties. Cultured glia, grown under the right conditions, also achieve characteristics similar to their *in vivo* counterparts. For example, oligodendrocyte

precursor cells can proliferate in culture and express normal oligodendrocyte markers. These precursor cells can also be differentiated into oligodendrocytes that produce two-dimensional “myelin sheets” on the surface of coverslips, allowing scientists to study the role of various proteins in forming a normal myelin sheath.

The main limitation of dissociated cell culture models is that organization and network connectivity of the *in vivo* environment is lost. Another limitation of dissociated cultures is the small quantity of cells relative to immortalized cell lines, which can make it difficult to perform biochemical experiments (Chapter 14) that require a high volume of starting material (millions of cells). Furthermore, most primary cell cultures are not homogeneous. Neuronal cultures are often mixtures of both glia and neurons that respond to different neurotransmitters, so identifying an individual population of cells can be difficult. To minimize heterogeneity, investigators usually attempt to dissect regions as precisely as possible to maximize the presence of the desired cell type in culture.

Various methods have been developed to purify specific cell populations. In a technique called **immunopanning**, a scientist coats the bottom of a plate with antibodies that recognize cell-surface markers on the outside of specific cell types. When heterogeneous populations of cells are added to the plate, the scientist can purify the cells of interest by allowing the cells to bind to the bottom of the plate and then wash off the undesired, unbound cells. In some cases, it is possible to achieve 95%–99% purity for a given cell type. This technique has been used to culture many different neural cell types, as well as oligodendrocyte precursors, astrocytes, and microglia.

STEM CELL CULTURES

A **stem cell** is a **pluripotent cell** with the capacity to generate any cell type (e.g., neuron, muscle, blood) and an unlimited ability to renew itself. These cells are a form of primary dissociated cells derived directly from a living organism. Yet they are also similar to immortalized cell lines in that they are theoretically able to propagate in culture indefinitely. Stem cells must be cultured in the presence of growth factors in specially formulated media to preserve their multipotency and capacity to self-renew.

Stem cells are classified by their source (e.g., embryonic stem [ES] cells, adult stem cells, induced pluripotent stem (iPS) cells), as well as the tissue they typically generate (e.g., neural stem cells [NSCs], hematopoietic stem cells, skin stem cells). Stem cells defined by the tissue they generate are **multipotent**—able to give rise to all types of cells found in the tissue and able to continuously self-renew. For example, NSCs can give rise to all three neural lineages—neurons, astrocytes, and oligodendrocytes—as well as additional NSCs. Stem cells are distinct from **progenitor cells**, which have a more limited capacity for self-renewal and may be **unipotent**, giving rise to a single cell type.

Stem cell culture has a great utility in neuroscience. ES cells are utilized for their ability to differentiate into specific neuronal subtypes. NSCs can be cultured to study the basic biology of development and aging. iPS cells can be used to generate clones of cells from patients with neurological disorders to characterize cellular and molecular changes in diseased neurons. While an extensive history and theory of stem cells already exists, basic stem cell biology is still a rapidly evolving field. Below, we describe specific categories of stem cells and their utility to neuroscience research.

Embryonic Stem Cells

Pluripotent **embryonic stem (ES) cells** can give rise to all tissues in an organism. ES cells are derived from the inner cell mass of a blastocyst, a structure that forms from a recently fertilized egg and develops into an embryo. Culturing ES cells essentially traps the cells in a pluripotent state by growing them in the presence of factors that prevent the cells from differentiating. Investigators can reconstruct the environment of ES cells in a dish so that culture conditions contain the specific molecules and mitogens that specify the formation of a neuron in the developing embryo.

Using special culture conditions, ES cells are first induced to become general **neural progenitors**, precursors that are committed to a neural fate. Once they have become neural progenitors, the specific molecules known to act during normal neuronal development are added to the culture to direct differentiation of specific neuron subtypes. For example, to make motor neurons, cells are exposed to high concentrations of sonic hedgehog (SHH) and retinoic acid, morphogens that pattern spinal cord motor neurons during development. To make midbrain dopaminergic neurons, cells are exposed to SHH and fibroblast growth factor 8, morphogens that specify dopaminergic fate during development.

Investigators monitor the progression of differentiation by examining the culture for activation of transcription factors relevant to development and the appearance of markers known to promote differentiation of specific neuronal subtypes. For example, when trying to induce midbrain dopaminergic neurons, scientists look for midbrain-specific transcription factors (Pitx3, En1, Lmx1b, Nurr1) and later for markers of mature dopamine neurons (TH, DAT, Girk2).

This process—from the initial establishment of the ES cells to their differentiation into specific neuronal subtypes—can take months. Also, despite attempts to create pure populations of a specific neuronal subtype, differentiated ES cell populations are still heterogeneous, with many other cells present in the culture. However, the ability to differentiate ES cells into specific neuronal cell types is useful for studying these cell types in an isolated, in vitro environment.

Neural Stem Cells

Multipotent **neural stem cells (NSCs)** exhibit self-renewal properties and the ability to differentiate into all neural subtypes. They can be extracted from regions of embryonic or adult brains, where they normally divide and give rise to neurons or glia. The embryonic brain contains many NSCs, while adult brains have far fewer NSCs. Furthermore, the ability of adult NSCs to produce neurons decreases with age.

Investigators typically culture adult NSCs from the subgranular zone of the dentate gyrus in the hippocampus and the subventricular zone in the wall of the lateral ventricles. After dissection, the region containing NSCs is dissociated and can be moderately purified through centrifugation. NSCs can also be grown and passaged in special media that promote NSC survival. By removing growth factors from the media, NSCs can be differentiated into neurons, astrocytes, and oligodendrocytes.

It is difficult to experimentally distinguish neural progenitors from NSCs, and scientists often confuse the terms. Both types of cells proliferate and express common sets of molecular markers (e.g., nestin, P_x6). Therefore, functional assays must also be used to determine the identity of cultured cells. Two defining features of NSC cultures that can be tested to confirm cellular identity are (1) the ability to give rise to all three types of the neural lineage (multipotency) and (2) the ability to propagate more NSCs (self-renewal).

Self-renewal is often examined in mice by performing **primary neurosphere** and **secondary neurosphere assays**. When NSCs are cultured in non-adherent conditions in the presence of the growth factors EGF and bFGF, they give rise to **neurospheres**, balls of dividing cells ([Fig. 13.4](#)). However,




FIGURE 13.4 **Neurospheres.** Dividing neural progenitor and stem cells form spheres in culture. *Courtesy of Dr. Victoria Rafalski.*

neurospheres form as a general result of dividing neural precursors, including progenitor cells that do not continuously self-renew. Therefore, while this assay is frequently used, it does not definitively distinguish between stem and progenitor cells. The secondary neurosphere assay involves culturing cells from established neurospheres to determine if the cells generated by the first neurosphere are able to continue proliferating, an indicator of self-renewal. These assays can be difficult to quantify, as cells and spheres can each fuse, so a sphere may not necessarily form from a single NSC. Also, neurospheres are species specific; rat and human NSCs do not form neurospheres as often as mouse NSCs.

Induced Pluripotent Stem Cells

Induced pluripotent stem (iPS) cells are differentiated cells that scientists have reverted back to a stem cell state. A variety of strategies have been employed to carry out this feat. One common strategy is to use viruses to introduce specific transcription factors into differentiated somatic cells, usually fibroblasts. The presence of these transcription factors seems to allow these cells to exhibit stem cell behaviors. Other strategies employ a combination of viral delivery and chemical manipulations to coax somatic cells into a pluripotent state.

Once iPS cells are generated, they are examined for characteristics of “stemness” by comparing them to ES cells in various ways. Scientists typically examine transcription factor profiles, measuring the presence of markers such as Oct4, Nanog, Sox2, AP, SSEA4, and TRA-1-80, using immunohistochemistry ([Chapter 6](#)) and reverse transcription polymerase chain reaction techniques ([Chapter 10](#)). The potential for these iPS cells to differentiate is measured by driving cells into specific lineages using established protocols.

In neuroscience, the ability to reprogram fully differentiated cells into an ES-like state opens up the possibility of modeling human disease in vitro, directly from patient-derived cells. These cells can reveal the cellular and molecular pathogenesis of diseases such as Parkinson’s and autism spectrum disorders. iPS cells may also be used for therapeutic purposes, for example, generating new neural cell types for transplantation from a patient’s own somatic cells. In some cases, it is even possible to bypass the iPS cell stage to directly reprogram somatic cells into functional neurons.

Brain Organoids

Brain organoids (also referred to as **cerebral organoids**) are derived from pluripotent stem cells and self-organize into 3D structures that can mimic the architecture and function of various brain regions, such as the cerebral cortex ([Fig. 13.5](#)). To generate brain organoids, researchers first transfer stem cell colonies onto low attachment plates to encourage the floating colonies to fold into spherical structures called **spheroids**. The spheroids are then pushed




FIGURE 13.5 Brain organoids systems: (A) Assembloids combine populations of organoids, each principally composed of a different cell type. (B) Organoids can be composed of many cell types to model specific parts of the brain (such as the cerebral cortex). (C) Spheroids with oligodendrocytes can be used to model *in vivo* neurodevelopment, especially the myelination of developing axons. (D) Because microglia do not develop in the brain itself, they can be harvested from other cell culture models and implanted into organoid cultures. (E) To study organoids in the context of the brain vasculature, organoids grown in culture can be grafted into rodent brains to promote vascularization. *Reprinted by permission from Oliveira, B. et al. (2019). Modeling cell-cell interactions in the brain using cerebral organoids. Brain Research: Elsevier.*

toward a neural fate via the addition of neutralizing growth factors. Depending on what part of the central nervous system the researcher wants to generate, different factors are applied to the spheroids to promote a particular regional fate. Protocols for the generation of organoids can take up to 70 days *in vitro*. In fact, researchers may have to wait up to 180 days for the development of astrocytes. Depending on the question, scientists may find it useful to create organoids that include glial cells such as oligodendrocytes, astrocytes, or microglia. Alternatively, many groups have been using assembloids, which are created from the combination of different types of patterned organoids. These advances in organoid complexity allow researchers to study neurons and glia in ways that were previously impossible. Organoids also allow scientists the ability to study human diseases in a more complex *in vitro* model than immortalized cell lines or primary cell culture allows.

MANIPULATING CELLS IN CULTURE

One of the major advantages of studying the nervous system *in vitro* is the ability to control the cellular environment. In addition to other methods of neuromodulation ([Chapter 8](#)), scientists can manipulate neurons in culture using several methods.

- **Transfection and Infection.** **Transfection** and **infection** are the processes of delivering DNA to cells using nonviral and viral methods, respectively ([Chapter 11](#)). One of the major benefits of cell culture systems, especially dissociated cultures, is that they are very amenable to both DNA delivery methods. Therefore, it is usually easy to introduce recombinant DNA molecules into cultured neurons to change gene/protein expression and alter cellular physiology.

- **Coculture Systems.** **Coculturing** allows a variety of cell types to be cultured together to examine the effect of one culture system on another (Fig. 13.6). This procedure is useful when examining the effect of one type of tissue on another, one region of the brain on another, or how a particular secreted molecule leads to changes in neural development or physiology. For example, cocultures of different regions of spinal cord explants initially revealed differing effects on the ability to attract or repel neurite outgrowth. Biochemical purification from explants in coculture experiments led to the identification of specific molecules that could then be introduced into immortalized cell lines to express and secrete these molecular guidance cues. Transfected cell lines could then be cocultured with spinal cord explant cultures to examine the neurite outgrowth response of the spinal cord neurons to specific guidance cues. Such approaches have deciphered a large variety of chemoattractants and chemorepellants that have since been validated *in vivo*.
- **Pharmacology.** Because of the unparalleled access to the cell culture media, it is relatively easy to add pharmacological agents to culture media to affect proteins, ion channels, and receptors of interest. Unlike *in vivo* conditions, it is possible to quickly wash out pharmacological reagents from the media to study the time course of removing the compound. Cultured cells are also very amenable to chemogenetic experiments in




FIGURE 13.6 The use of cocultures to study guidance cues. In this example, midbrain explants are placed in culture medium. (A) COS cells (top) transfected with a control plasmid do not influence the direction of axon outgrowth from the explants. (B) COS cells transfected with netrin-1 attract axons from the explants. *Courtesy of Dr. Jie Li and Dr. Mary Hynes.*

which a synthetic drug causes a change in activity by binding to a synthetic receptor. See [Chapter 8](#) for a detailed description of pharmacological and chemogenetic techniques.

- *Optogenetics.* Optogenetic techniques ([Chapter 8](#)), in which neurons are transduced with a light-activated transmembrane protein to cause changes in neural activity, are relatively easy to perform in culture conditions. Light can be delivered from a nearby fiber optic cable or can even be delivered directly through the microscope objective near the cell culture plate. Optogenetic neuromodulation of cells in culture is often combined with electrophysiological recordings to study the effects of optogenetic actuators in different cell types or to study circuit connections in cultured brain slices.
- *Antibody Interference.* Antibodies are molecules that bind with high affinity to antigens, usually parts of a protein. In nature, these antibodies function as part of the immune system. In the lab, antibodies can be used in many different assays to study protein expression and binding partners ([Chapter 14](#)). If the antigen of an antibody is accessible in living cells, it is also sometimes possible to block protein function by applying antibodies against the protein of interest. These antibodies bind the protein at a site critical for its normal function, thus inhibiting function in vitro or in vivo.

CONCLUSION

The ability to grow cells and tissue outside of a living organism allows scientists much greater access and control over the cellular environment compared with in vivo conditions. Therefore, scientists can investigate the cellular and molecular mechanisms of nervous system function with great detail. The tools and techniques of in vitro culture techniques can provide many informative results, but ultimately this knowledge must be tested in the intact nervous system to fully understand how a biological process works in a living organism.

SUGGESTED READINGS AND REFERENCES

Books

- Barres, B.A., Stevens, B. (Eds.), 2013. *Purifying and Culturing Neural Cells: A Laboratory Manual*. Cold Spring Harbor Laboratory Press, New York.
- Butler, M., 2004. *Animal Cell Culture and Technology*, second ed. BIOS Scientific Publishers, London.
- Freshney, R.I., Capes-Davis, A., 2021. *Freshney's Culture of Animal Cells: A Manual of Basic Techniques and Specialized Applications*, fifth ed. Wiley-Liss, Hoboken, NJ.
- Skloot, R., 2011. *The Immortal Life of Henrietta Lacks*. Broadway Books, New York.

Review Articles

- Amin, N.D., Pasca, S.P., 2019. Building models of brain disorders with three-dimensional organoids. *Neuron* 100, 389–405.
- Blau, H.M., Brazelton, T.R., Weimann, J.M., 2001. The evolving concept of a stem cell: entity or function? *Cell* 105, 829–841.
- Jaenisch, R., Young, R., 2008. Stem cells, the molecular circuitry of pluripotency and nuclear reprogramming. *Cell* 132, 567–582.
- Marton, R.M., Pasca, S.P., 2020. Organoid and assembloid technologies for investigating cellular crosstalk in human brain development and disease. *Trends Cell Biol.* 30, 133–143.
- Oliveira, B., et al., 2019. Modeling cell-cell interactions in the brain using cerebral organoids. *Brain Res.* 1724, 146458.

Primary Research Articles—Interesting Examples from the Literature

- Barberi, T., Klivenyi, P., Calingasan, N.Y., Lee, H., Kawamata, H., Loonam, K., Perrier, A.L., Bruses, J., Rubio, M.E., Topf, N., Tabar, V., Harrison, N.L., Beal, M.F., Moore, M.A., Studer, L., 2003. Neural subtype specification of fertilization and nuclear transfer embryonic stem cells and application in parkinsonian mice. *Nat. Biotechnol.* 21, 1200–1207.
- Chiaradia, I., Lancaster, M.A., 2020. Brain organoids for the study of human neurobiology at the interface of in vitro and in vivo. *Nat. Neurosci.* 23, 1496–1508.
- Dotti, C.G., Banker, G.A., 1987. Experimentally induced alteration in the polarity of developing neurons. *Nature* 330, 254–256.
- Dugas, J.C., Mandemakers, W., Rogers, M., Ibrahim, A., Daneman, R., Barres, B.A., 2008. A novel purification method for CNS projection neurons leads to the identification of brain vascular cells as a source of trophic support for corticospinal motor neurons. *J. Neurosci.* 28, 8294–8305.
- Foo, L.C., Allen, N.J., Bushong, E.A., Ventura, P.B., Chung, W.S., Zhou, L., Cahoy, J.D., Daneman, R., Zong, H., Ellisman, M.H., Barres, B.A., 2011. Development of a method for the purification and culture of rodent astrocytes. *Neuron* 71, 799–811.
- Guttenplan, K.A., Weigel, M.K., Prakash, P., Wijewardhane, P.R., Hasel, P., Rufen-Blanchette, U., Munch, A.E., Blum, J.A., Jine, J., Neal, M.K., Bruce, K.D., Gitler, A.D., Chopra, G., Liddelow, S.A., Barres, B.A., 2021. Neurotoxic reactive astrocytes induce cell death via saturated lipids. *Nature* 599, 102–107.
- Pan, P.Y., Tian, J.H., Sheng, Z.H., 2009. Snapin facilitates the synchronization of synaptic vesicle fusion. *Neuron* 61, 412–424.
- Paşca, A., Sloan, S., Clarke, L., et al., 2015. Functional cortical neurons and astrocytes from human pluripotent stem cells in 3D culture. *Nat. Methods* 12, 671–678.
- Rafalski, V.A., Ho, P.P., Brett, J.O., Ucar, D., Dugas, J.C., Pollina, E.A., et al., 2013. Expansion of oligodendrocyte progenitor cells following SIRT1 inactivation in the adult brain. *Nat. Cell Biol.* 15, 614–624.
- Raineteau, O., Rietschin, L., Gradwohl, G., Guillemot, F., Gähwiler, B.H., 2004. Neurogenesis in hippocampal slice cultures. *Mol. Cell Neurosci.* 26, 241–250.
- Soldner, F., Hockemeyer, D., Beard, C., Gao, Q., Bell, G.W., Cook, E.G., Hargus, G., Blak, A., Cooper, O., Mitalipova, M., Isacson, O., Jaenisch, R., 2009. Parkinson's disease patient-derived induced pluripotent stem cells free of viral reprogramming factors. *Cell* 136, 964–977.
- Takahashi, K., Yamanaka, S., 2006. Induction of pluripotent stem cells from mouse embryonic and adult fibroblast cultures by defined factors. *Cell* 126, 663–676.

- Tessier-Lavigne, M., Placzek, M., Lumsden, A.G., Dodd, J., Jessell, T.M., 1988. Chemotropic guidance of developing axons in the mammalian central nervous system. *Nature* 336, 775–778.
- Vierbuchen, T., Ostermeier, A., Pang, Z.P., Kokubu, Y., Sudhof, T.C., Wernig, M., 2010. Direct conversion of fibroblasts to functional neurons by defined factors. *Nature* 463, 1035–1041.
- Watanabe, K., Kamiya, D., Nishiyama, A., Katayama, T., Nozaki, S., Kawasaki, H., Watanabe, Y., Mizuseki, K., Sasai, Y., 2005. Directed differentiation of telencephalic precursors from embryonic stem cells. *Nat. Neurosci.* 8, 288–296.

Protocols

- Current Protocols in Neuroscience. Chapter 3: Cellular and Developmental Neuroscience, 2006. John Wiley & Sons, Inc., Hoboken, NJ.
- Loring, J.F., Wesselschmidt, R.L., Schwartz, P.H. (Eds.), 2007. *Human Stem Cell Manual: A Laboratory Guide*. Academic Press, Amsterdam.
- Poindron, P., Piguet, P., Förster, E. (Eds.), 2005. *New Methods for Culturing Cells from Nervous Tissues*. Karger, Basel.
- Zigova, T., Sanberg, P.R., Sanchez-Ramos, J.R. (Eds.), 2002. *Neural Stem Cells: Methods and Protocols*. Humana Press, Totowa, NJ.

Websites

- American Type Culture Collection: <http://www.atcc.org>.
- European Collection of Cell Cultures: <http://www.hpacultures.org.uk/aboutus/ecacc.jsp>.
- Life Technologies Protocols: <http://www.lifetechnologies.com/us/en/home/references/protocols.html>.